

TouchNTuff® 83-300 83-500



Single Use



Polyisoprene



Splash



Polyisoprene Gloves. Comfortable Latex-Free gloves

DESCRIPTION

- Made of 100% polyisoprene, the TouchNTuff® 83-300 contains no natural rubber latex, minimising the risk of Type I allergy.
- The stretchability of the ultra-soft polyisoprene film and a new ergonomic former provide second-skin fit and greater comfort. Users can wear the glove longer before developing hand fatigue.
- The synthetic polyisoprene delivers all the tactile performance of natural rubber latex.
- An anti-slip finish provides superior grip when handling wet items and helps prevent slippage.
- The enhanced polymer coating inside the glove allows easy donning with dry or damp hands. Design and thickness allow for easy and comfortable double-gloving.

PRODUCT DETAILS

Coating material	Polyisoprene
Grip design	Smooth Finish
Cuff style	Rolled beaded
Thickness (mm)	83-300: 0,21 83-500: 0,21
Colour	White
AQL	1.0
EN Sizes	5.5,6,6.5,7,7.5,8,8.5,9
Length (mm)	305-310

PACKAGING

83-300: 20 right and left hand gloves packed in separate compartments in a polybag; 1 polybag per vacuum sealed polybag; 2 polybags per vacuum sealed master polybag, 5 master polybags per case liner; 200 pairs per carton. 83-500: Pair pack: The glove is packaged one pair per poly wallet, sealed in a poly envelope; 10 pairs per sealed polybag; 5 sealed polybags per master bag; 4 master bags per case liner; 200 pairs per carton. All packaging is printed with cleanroom-compatible IPA-resistant ink. The polybags have an easy-tear opening. The gloves are triple-bagged for superior contamination control.

IDEAL APPLICATIONS

- Analytical measurements
- Blending of liquids
- Filtration processes
- Handling bags containing 'active' & 'inactive' materials
- Pelletizing and tablet manufacturing
- Tablet press preparation and operation
- Vaccines injectables operations
- Wet chemical, particle and analytical testing
- 'V' Blender to mix solids and colorants

PRIMARY INDUSTRIES



STANDARDS AND CERTIFICATIONS



CATEGORY III

FEATURED TECHNOLOGY

EXTRA FEATURES

Print date 19-06-2019

Ansell Healthcare Europe N.V. (European Head Office)

Riverside Business Park, Block J

Boulevard International 55, 1070 Brussels, Belgium

Tel. +32 (0) 2 528 74 00 • Fax +32 (0) 2 528 74 01 • Fax Customer Service +32 (0) 2 528 74 03

<http://www.ansell.eu> • E-mail info@ansell.eu

